

ULTRA™

Ultra power.
Ultra utility.

CYNOSURE®
LUTRONIC®





The most versatile non-ablative *1927nm Fractional Laser*

The highest-powered (20W) 1927nm fractional laser effectively rebuilds glowing and healthy-looking skin with minimal downtime.

With the ability to lighten sun damage and increase superficial collagen, Ultra enables easy customization so all patient skin types can be treated all year round.





Patient and practice *Benefits*



Versatility to treat all skin types:
From pre-juvenating and maintaining
a youthful complexion to more
aggressive revitalization.



Treatment add-ons available with topicals,
such as tranexamic acid for lighter
skin tone, and specialized cream masks
for enhanced post-treatment experience.



Smart IntelliTrak™ Rolling System
for rapid and consistent energy delivery.



Controlled Chaos Technology (CCT)
ensures intelligent, rapid placement
of microbeams to reduce risk of bulk
heating and patient discomfort.



Measured energy delivery for enhanced
clinical confidence.



New, complementary treatments
with Ultra Skin Solutions can provide
exceptional ROI.

Real patients. Genuine concerns. *Life-changing results* /

Courtesy of Kate Holcomb, MD



Before and after 3 treatments

Courtesy of James Y. Wang, MD



Before and after 2 treatments

Courtesy of Lemmon Avenue



Before and after 2 treatments

Courtesy of Kate Holcomb, MD



Before and after 2 treatments



ULTRA+

SKIN SOLUTIONS

Immediate Results, Lasting Benefits.

Ultra can perform complementary treatments with various types of devices, such as full beam pigment lasers and RF microneedling.

With RF Microneedling:

The treatment targets multiple depths for complementary multi-level results. Microneedling rebuilds the structure and foundation of the skin, stimulating collagen development, while the Ultra treatment delivers fast and effective non-ablative revitalization of superficial skin layers.

With Full Beam Pigment Lasers:

Ultra can uniquely complement full beam pigment laser treatments to address common patient concerns effectively and in fewer sessions. Together, they deliver immediate smoothing and brightening, as well as lasting collagen building.





Topicals and treatment *enhancements*

Ultra laser treatments create tiny dehydrated “microchannels” in the skin which increase absorption of specially curated topicals designed to address patient concerns.

LaseMD ampoules (vitamin A, vitamin C, tranexamic acid, and resveratrol) are designed to enhance skin revitalization. The ampoules use pure ingredients with no perfumes, dyes, or preservatives, and they are specifically formulated nano particles to help with absorption. Use of a bio cellulose cream mask can optimize patient comfort and results, and patients love the moisturizing, soothing and cooling sensation post treatment. Ampoules and masks make for an exceptional patient experience.





We're your trusted *partner on this journey/*

Technology

- Controlled Chaos Technology
- Forward thinking science
- World's highest powered 1927nm laser

Product

- Significant clinical evidence
- Highest quality technological improvements

Practice Development

- Industry expert consultants
focused on growing your business

Clinical

- Comprehensive virtual training
- Personalized in-office training for beginner
and advanced teams

Service

- Experienced technical support team 24/7
- Expert field service team

Amplified Marketing Practice Support

- Dedicated patient support team
- Comprehensive turnkey marketing



ULTRA™

Visit/
www.lutronic.com

Follow us:

 @cynosurelasers

 Cynosure, LLC

 @lutronic_us

 Lutronic

Intended for distribution in US only.

The Ultra aesthetic workstation is indicated for the following:
Intended for use in dermatological procedures requiring the coagulation of soft tissue;
Treatment of actinic keratosis; Treatment of benign pigmented lesions; Lentigos (age spots); Solar lentigos (sun spots); Ephelides (freckles).

Cynosure, LLC and Lutronic Aesthetic, Inc. do not diagnose or treat medical conditions, prescribe treatment protocols, provide medical advice, or otherwise practice medicine. Patient should consult with a qualified healthcare professional for diagnosis and proper treatment protocols/parameters to determine if treatment with any product is right for them. Qualified healthcare professionals are responsible for diagnosis and treatment protocols using products.

©2024 CYNOSURE LUTRONIC, INC. All Rights Reserved. LaseMD Ultra, Controlled Chaos Technology are trademarks of Lutronic Aesthetics. PRD-1434

921-7068-001 R01



Attributes & Specifications

Attributes	Specifications
Laser Type	1927nm Thulium
Platform Power	20W
Laser Delivery System	Optical Fiber with Handpiece
Pulse Energy	1 to 20mJ (1mJ increments)
Pulse Rate	Roller, Square, Comb Tips: 143-323Hz C5 Tip: 182-632Hz
Pulse Duration	Roller, Square, Comb Tips: 0.1-4ms C5 Tip: 0.083 – 4ms
Controlled Chaos Technology™	Exclusively available to reduce risk of bulk heating and PIH
Ultra Tips Included	Roller Tip – 200µm diameter microbeams C1 Tip – 200µm diameter microbeams Comb Tip – 200µm diameter microbeams for KeraLase C5 Tip – 350µm diameter microbeams for discreet pigment
Dimensions and Weight	15.4" W x 24.1" D x 68.7" H & 72.8lbs (33kg)
Electrical Rating	150VA (AC 100-240V) 50/60Hz. Fuse 250V/6.3A

